

Test Definition Report

09/30/2015 09:18:34

✓ **Enabled** Test code MRX143 Test name D-dimer MRX143 Test number 500

General information

Test code Test name
LIS number Test number

Lock manual ordering

Assay

Assay release number

Test modification number

User revision comment

IL revision comment

Manual validation required

Analytical cycle definition

Acquisition - wavelength

Primary wavelength nm

Acquisition - time

Standard time s

Delay time s

Load cycle

Number of replicates

Material Name	Type	Load Volume	Incubation Range
Sample		15	20.000 - 55.000
D-dimer RB MRX143	Intermediate reagent	95	20.000 - 55.000
D-dimer Latex MRX143	Start reagent	55	

Total cuvette volume μ L

Test Definition Report

09/30/2015 09:18:34

Enabled Test code **MRX143** Test name **D-dimer MRX143** Test number **500**

Material/Sample Load Cycle Definition

Material

Aspiration information

Head volume	<input type="text" value="0"/> μL	Material/Sample volume	<input type="text" value="15"/> μL
Airgap	<input type="text" value="15"/> μL	Transport airgap	<input type="text" value="15"/> μL
Head volume airgap	<input type="text" value="0"/> μL	Total volume	<input type="text" value="45"/> μL

Dispensation information

Enhanced dispense
Enable mix

Incubation range s

Before Aspiration Rinse / Clean & Rinse

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

Lock rinse / clean & rinse

Rinse

Time s
Enable agitation

Clean & rinse

Material

Aspiration information

Head volume	<input type="text" value="0"/> μL	Material/Sample volume	<input type="text" value="95"/> μL
Airgap	<input type="text" value="15"/> μL	Transport airgap	<input type="text" value="15"/> μL
Head volume airgap	<input type="text" value="0"/> μL	Total volume	<input type="text" value="125"/> μL

Dispensation information

Enhanced dispense
Enable mix

Incubation range s

Before Aspiration Rinse / Clean & Rinse

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

Rinse

Time s Frequency
Enable agitation

Clean & rinse

Material

Aspiration information

Head volume	<input type="text" value="0"/> μL	Material/Sample volume	<input type="text" value="55"/> μL
Airgap	<input type="text" value="15"/> μL	Transport airgap	<input type="text" value="15"/> μL
Head volume airgap	<input type="text" value="0"/> μL	Total volume	<input type="text" value="85"/> μL

Dispensation information

Enhanced dispense
Enable mix

Before Aspiration Rinse / Clean & Rinse

Test Definition Report

09/30/2015 09:18:34

Enabled Test code MRX143 Test name D-dimer MRX143 Test number 500

Material/Sample Load Cycle Definition

Material

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

Rinse

Time s Frequency

Enable agitation

Clean & rinse

Enable sample pre-dilution

Alternative pre-dilution

Use extended acquisition time

Diluted LIS number

Test code

Dilution definition

Number of parts sample	<input type="text" value="1"/>	
Number of parts diluent	<input type="text" value="14"/>	
Minimum aspiratable volume	<input type="text" value="15"/> μL	
Volumes (μL)		
Sample	Mixture	Diluent
15		210

Enable mix

Mix volume	<input type="text" value="50"/> %
Number of cycles	<input type="text" value="2"/>
Rinse after mix	<input type="text" value="1"/> s

Sample/Mixture

Volumes

Head volume	<input type="text" value="0"/> μL
Airgap	<input type="text" value="15"/> μL
Transport airgap	<input type="text" value="15"/> μL
Head volume airgap	<input type="text" value="0"/> μL

Before Aspiration Rinse / Clean & Rinse

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

Rinse

Time s

Enable agitation

Clean & rinse

Test Definition Report

09/30/2015 09:18:34

✓ **Enabled** Test code **MRX143** Test name **D-dimer MRX143** Test number **500**

Diluent

Volumes

Material	Factor	Diluent
Head volume	0	µL
Airgap	15	µL
Transport airgap	15	µL
Head volume airgap	0	µL

Before Aspiration Rinse / Clean & Rinse

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

✓ Rinse

Time	<input type="text" value="1"/> s	Frequency	<input type="text" value="Between changes in material only"/>
Enable agitation	<input type="text" value="1"/>		

Clean & rinse

✓ DR Parameters - Primary wavelength - Wavelength definition

Primary wavelength: Enable warnings Enable errors

✓ Enable normalization method

Smoothing

- ✓ Enable smooth 1 Number of points
- ✓ Enable smooth 2 Number of points
- ✓ Enable Smooth 3 Number of points

Enable raw data checks

✓ Enable normalized data checks

Baseline Check

✓ Enable baseline check Number of seconds

✓ Enable primary algorithm

Linear kinetic

✓ Slope of curve Start time s Time s

✓ Enable linear regression SD error check SD error limit

✓ Enable secondary algorithm

Linear kinetic

✓ Delta of curve Method
Start time s
End time s

✓ Enable negative calculation

✓ Enable initial slope check Time s
Maximum error limit Maximum warning limit

Result units

Unit selection Millabsorbance/min mAbs/min

Unit	<input type="text" value="Millabsorbance/min"/>	Unit type	<input type="text" value="Measured"/>
Label	<input type="text" value="mAbs/min"/>	Decimals	<input type="text" value="1"/> IL Unit
User defined label			

✓ Enable normal range Range

Test Definition Report

09/30/2015 09:18:34

✓ **Enabled** Test code **MRX143** Test name **D-dimer MRX143** Test number **500**

Unit selection **Milliabsorbance/min** **mAbs/min**

Allow negative results (QC jobs excluded)

Unit selection **ng/mL** **ng/mL**

Unit	ng/mL	Unit type	Calibrated
Label	ng/mL	Decimals	0 IL Unit
User defined label			

✓ Enable test range

Range	50 - 4000	Extended range	-999999 - 999999
Alternate range	0 - 0		

Allow negative results (QC jobs excluded)

Primary unit cfg.

Primary unit **ng/mL**

Patient result selection

Unit 1	ng/mL
Unit 2	
Unit 3	
Unit 4	

Measured result curve display settings

<input checked="" type="checkbox"/> Enable auto scale	
<input type="checkbox"/> Raw	<input checked="" type="checkbox"/> Normalized
Y-Axis minimum	Y-Axis maximum
0.000	3000.000 mAbs

Calibration

✓ Enable calibration

General

Number of replicates	4
Calibration unit	ng/mL
Calibration mode	<input checked="" type="checkbox"/> Automatic dilutions process <input type="checkbox"/> Prediluted calibrators

Test not feasible when no validated calibration

Force mode with reduced clean

Enable import calibration

Frequency

Enable frequency

DR Parameters

✓ Enable extrapolation	Lower limit 1.000 Upper limit 3600.000
------------------------	--

✓ Enable outlier removal

Math Model

Math model	Spline	First concentration	1
X Axis transformation	No transformation	Last concentration	5
Y Axis transformation	No transformation		

✓ Automatic dilution

<input checked="" type="checkbox"/> Define points by % dilution	<input type="checkbox"/> Define points by target concentration
---	--

Dilution process

<input type="checkbox"/> Direct	<input checked="" type="checkbox"/> Serial
<input type="checkbox"/> Single dilution	<input checked="" type="checkbox"/> Batch

Calibrator target value **2755.000** ng/mL

Test Definition Report

09/30/2015 09:18:34

✓ **Enabled** Test code **MRX143** Test name **D-dimer MRX143** Test number **500**

DR Checks Configuration
Raw Data Check Def Norm. Data Check Def.

Conc.	Target value	Min. Replic	% CV Max.	%Max. mean	%Min. mean	Volumes (µL)			Predilution volumes (µL)			First Point	Last Point	Signal	Curve Sequence	Enable spike removal	Baseline By Moving SD	Norm. Delta Upp. Lim.	Norm. Delta Low. Lim.	Calibrated Result Delta	Baseline	Endpoint	Multiple Threshold
						Calibrator	Mixture	Diluent	Sample	Mixture	Diluent												
100.00	2755	1	5.000			243																	
50.00	1378	1	5.000	100.000	0.000		118	118															
25.00	689	1	5.000	100.000	0.000		110	110															
12.50	344	1	5.000	100.000	0.000		94	94															
6.25	172	1	5.000	100.000	0.000		63	63															

Minimum aspiratable volume

✓ Enable mix

Mix volume	<input type="text" value="50"/>	%
Number of cycles	<input type="text" value="2"/>	
Rinse after mix	<input type="text" value="1"/>	s

Calibrator

Volumes

Material	<input type="text" value="DD Calibrator3200"/>
Head volume	<input type="text" value="0"/> µL
Airgap	<input type="text" value="15"/> µL
Transport airgap	<input type="text" value="15"/> µL
Head volume airgap	<input type="text" value="0"/> µL

Before Aspiration Rinse / Clean & Rinse

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

✓ Rinse

Time	<input type="text" value="1"/>
Enable agitation	<input type="checkbox"/>

Clean & rinse

Diluent

Volumes

Material	<input type="text" value="Factor Diluent"/>
Head volume	<input type="text" value="0"/> µL
Airgap	<input type="text" value="15"/> µL
Transport airgap	<input type="text" value="15"/> µL
Head volume airgap	<input type="text" value="0"/> µL

Before Aspiration Rinse / Clean & Rinse

Enable rinse / clean & rinse before aspiration

After Dispensation Rinse / Clean & Rinse

✓ Rinse

Time	<input type="text" value="1"/>	Frequency	<input type="text" value="Between changes in material only"/>
Enable agitation	<input type="checkbox"/>		

Clean & rinse

Test Definition Report

09/30/2015 09:18:34

✓ **Enabled** **Test code** MRX143 **Test name** D-dimer MRX143 **Test number** 500

Enable parallelism

Enable NPP

✓ **Enable rerun**

Result Ranges

Result unit to test

ng/mL

Rule

Same Test

Extended Test

Alternative pre-dilution

Above Test Range

✓

Below Test Range

✓